





Australian Plant Phenomics Facility


The High Resolution Plant Phenomics Centre



Data acquisition and management

Dr Peter Ansell

Overview



Data acquisition

- **Experimental layouts determined before experiment begin**
- **Layout uploaded to machine databases and PODD**
- **Images and other measurements collected regularly using PlantScan, TrayScan, CabScan, Cropatron, etc.**
- **Images stored on disks specific to each platform, and metadata stored in SQL databases for each platform**

Data management and processing

- **Images are batched up and sent to the CSIRO Datastore**
- **CSIRO Datastore is backed up to a long term tape library**
- **Images are processed on the CSIRO Pearcey cluster and in the near future on the NCI node of the Nectar cloud**
- **Results are sent to the CSIRO Datastore and backed up on tape**

Data analysis

- **Results are aggregated based on the experimental layout that is found in PODB**
- **Categories include:**
 - **Genotype**
 - **Pot**
 - **Treatment**
 - **Replicate**
- **For each category, the average, maximum, minimum, standard deviation and count are computed for each day where data was acquired**

Discussion

- **Experimental layouts must be accurate for analysis to be useful**
- **CSIRO Datastore capacity makes it possible to permanently store our very large datasets**
- **Experimenting with new technologies**

Future work

- **Further visualisation of data and results**
- **Integration with genomics**
- **Expand the number of platforms that use the CSIRO Datastore for permanent archival**
- **Expand the number of platforms that use cluster processing**
- **Define data standards to enable other scientists to send us their images and experimental layouts for us to attempt to process them**